JAUNE	ALLIAGE DENTAIRE CONVENTIONNEL
GELB	DENTAL GÜSSELEGIERUNGEN
GIALLO	LEGA DENTALE CONVENZIONALE
YELLOW	DENTAL CASTING ALLOY
AMARILLO	ALEACION DENTAL CONVENCIONAL

% = Au: 51.0, Pd: 7.0, Ag: 27.0, Cu: 14.0, Zn, Ir : < 1%

13.4 g/cm²
870-920°C
290 HV1

ISO 22674

Fabriqué par / Hergestellt von / Fabbricato da /
Made by / Fabricado por

PX DENTAL SA

Champs-Montants 16a, 2074 Marin

Switzerland

T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29

www.pxdental.com

PXDENTAL SA

C H - 2 0 7 4 M A R I N

Mode d'emploi

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1. Waxing up, spruing, investing

- Follow the general prescriptions for casting alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- All regular or phosphate-bonded investments for precious metal alloys may be used. For the alloy PX Gold 51 we recommend to use the PX Extravest investment.

2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 680° C for a soaking period of 30-45 min.

3. Crucibles and casting buttons

- PX Gold 51 may be cast in a graphite, quartz or ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

4. Melting

- For PX Gold 51 casting temperature is 1070° C with a graphite crucible, 1120° C with a ceramic crucible.
- The use of a pinch of Borax is recommended.
- Bench cool the casting ring.

5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Finishing: use abrasive grinding points of your choice. Clean the frameworks in an ultrasonic unit or with a steam cleaner.

6. Polishing

- For best prepolishing results use soft brushes, felts and cotton wheels with a diamond paste. For final high polish use a polishing paste or "whiting chalk" mixed with water on soft brushes or cotton wheels.

7. Soldering

- Apply the solder PX SOLD 830, yellow, at 830° C using a propane/ oxygen torch and a pickling paste.
- Apply the solder PX SOLD 780, yellow, at 780° C using a Propane /oxygen torch and a pickling paste.

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Fabriqué par / Hergestellt von / Fabbricato da / Made by / Fabricado por

PX DENTAL SA

Champs-Montants 16 a, CH-2074 Marin
Switzerland

T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29

www.pxdental.com