BLANC	ALLIAGE DENTAIRE CONVENTIONNEL
WEISS	DENTAL GUSSLEGIERUNGEN
BIANCO	LEGA DENTALE CONVENZIONALE
WHITE	DENTAL CASTING ALLOY
BLANCO	ALEACION CONVENCIONAL

% = Au: 3.0, Pd: 23.0, Ag: 65.0, Cu: 6.5, In 2.0, Ir, Zn: < 1%

10.7 g/cm<sup>2</sup>  
990-1050°C  
180 HV1

ISO 22674

Fabriqué par / Hergestellt von / Fabbricato da /  
Made by / Fabricado por

### PX DENTAL SA

Champs-Montants 16a, 2074 Marin

Switzerland

T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29

[www.pxdental.com](http://www.pxdental.com)

**PXDENTAL SA**

C H - 2 0 7 4 M A R I N

## Mode d'emploi

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# Instructions for use

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## 1. Waxing up, spruing, investing

- Follow the general prescriptions for casting alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- All regular or phosphate-bonded investments for precious metal alloys may be used. For the alloy PX Ecocast G3 we recommend to use the PX Extravest investment.

## 2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 680° C for a soaking period of 30-45 min.

## 3. Crucibles and casting buttons

- PX Ecocast G3 may be cast in a graphite, quartz or ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

## 4. Melting

- For PX Ecocast G3 casting temperature is 1240°C with a ceramic crucible.
- The use of a pinch of Borax is recommended.
- Bench cool the casting ring.

## 5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Finishing: use abrasive grinding points of your choice. Clean the frameworks in an ultrasonic unit or with a steam cleaner.

## 6. Polishing

- For best prepolishing results use soft brushes, felts and cotton wheels with a diamond paste. For final high polish use a polishing paste or "whiting chalk" mixed with water on soft brushes or cotton wheels.

## 7. Soldering

- Apply the solder PX SOLD 830, yellow, at 830° C using a propane/ oxygen torch and a pickling paste.
- Apply the solder PX SOLD 780, yellow, at 780° C using a Propane /oxygen torch and a pickling paste.

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Fabriqué par / Hergestellt von / Fabbricato da / Made by / Fabricado por

**PX DENTAL SA**

Champs-Montants 16 a, CH-2074 Marin  
Switzerland

T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29

[www.pxdental.com](http://www.pxdental.com)